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
1. Louise M. Sternicki, Jim Nonomiya, Miaomiao Liu, et al. Native Mass Spectrometry for the Study of PROTAC GNE-987-Containing Ternary Complexes [J]. Chem Med Chem, 2021, 16, 1-6.
2. Lawrence B. Snyder, ArvinasInc, New Haven, CT. Discovery of ARV-110, a first in class androgen receptor degrading PROTAC® for the treatment of men with metastatic castration resistant prostate cancer. www.arvinas.com

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Quantitative Analysis of the molecular glue drugs Eragidomide in Plasma by SCIEX Triple Quad™ System

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1. Surka C, Jin L, Mbong N, et al. CC-90009, a novel cereblon E3 ligase modulator, targets acute myeloid leukemia blasts and leukemia stem cells[J]. Blood, 2021, 137(5): 661-677.
2. Hansen J D, Correa M, Alexander M, et al. CC-90009: a cereblon E3 ligase modulating drug that promotes selective degradation of GSPT1 for the treatment of acute myeloid leukemia[J]. Journal of Medicinal Chemistry, 2021, 64(4): 1835-1843.
3. Dong G, Ding Y, He S, Sheng C. Molecular Glues for Targeted Protein Degradation: From Serendipity to Rational Discovery. J Med Chem. 2021 Aug 12;64(15):10606-10620. doi: 10.1021/acs.jmedchem.1c00895. Epub 2021 Jul 28. PMID: 34319094.

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Improved LC-MRM Quantification Sensitivity for Powder for Inhalation from the Beta2 Agonist Family

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Key Words: Quantification, Indacaterol, Olodacaterol, Beta2 Agonist, Inhalation, MRM, Triple Quad

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1. <https://www.who.int/ncds/en>.
2. Wang C, et al. Prevalence and risk factors of chronic obstructive pulmonary disease in China (the China Pulmonary Health [CPH] study): a national cross-sectional study[J]. Lancet. 2018, 391(10131):1706-1717.
3. Global Initiative for Chronic Obstructive Lung Disease. Global strategy for the diagnosis, management and prevention of chronic obstructive pulmonary disease 2020 report[J]. <https://goldcopd.org/gold-reports/>
4. Vaidya S, et al. Pharmacokinetics of indacaterol, glycopyrronium and mometasone furoate following once-daily inhalation as a combination in healthy subjects. Pulm Pharmacol Ther 2020 10;64.

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Determination of Paxlovid on SCIEX Triple Quad™ LC-MS/MS

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Key words: SCIEX Triple Quad™ LC-MS/MS; Paxlovid

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Achieving low-pg/mL quantitation of leuprolide in human plasma

This technical note demonstrates a sensitive method to quantify leuprolide in human plasma on a SCIEX 7500 system. A lower limit of quantitation (LLOQ) of 2.5 pg/mL was determined using a 5-minute LC-MS/MS method (Figure 1).

Leuprolide acetate is a synthetic peptide GnRH agonist that inhibits the synthesis of estrogen and androgen hormones in women and men (Figure 1B). Long-acting depot formulations have been used in suppressive treatments in men with prostate cancer, women with endometriosis and uterine fibroids and pediatric patients with central precocious puberty.¹ Aside from long-acting depot formulations, leuprolide could be administered daily in microdose regimens.

Since different dosage regimens are required in hormone therapy, sensitive and selective assays for high-confidence detection and quantitative performance in biological matrices are needed to ensure the safety and efficacy of the promising therapeutic peptides.

Quantitation of peptides and proteins has been commonly performed using immunoassays, such as ELISA. However, immunoassays inherently face analytical challenges, such as poor selectivity, limited dynamic range and cross-reactivity.

LC-MS/MS methods are increasingly applied for the quantitation of therapeutics as they offer the most sensitive and selective platforms for bioanalysis.

The presented method demonstrates a sensitive method to quantify leuprolide in human plasma using an efficient solid-phase extraction (SPE) technique. Sensitivity at the low-pg/mL level was achieved across a linear dynamic range (LDR) spanning 3.6 orders of magnitude in human plasma using the SCIEX 7500 system. Front-end enhancements on the instrument facilitated improved overall ion generation, capture and transmission, therefore providing optimal quantitative sensitivity.

Key features of the quantitation of leuprolide using the SCIEX 7500 system

- Low-pg/mL level quantitation of a GnRH-releasing agonist: Achieve a 2.5 pg/mL LLOQ for leuprolide in human plasma on the SCIEX 7500 system
- Ideal analytical performance: Achieve accurate quantitative performance with %CV <10% at all


Figure 1. Representative extracted ion chromatograms (XICs) and structure of leuprolide. An LLOQ of 2.5 pg/mL (0.025 pg on column) was achieved for leuprolide in human plasma. No matrix interference was observed at the retention time of the analyte (A). Leuprolide is a synthetic nonapeptide used to suppress gonadal steroid synthesis (B).

concentration levels across a linear dynamic range (LDR) spanning 3.6 orders of magnitude

- Enhanced sensitivity unlocked: Improved front-end technology with the D Jet ion guide, OptiFlow Pro ion source and E Lens probe enhanced ion generation, capture and transmission, enabling users to reach desired quantitative sensitivity
- Streamlined data management: Data acquisition and processing are integrated into SCIEX OS software, a 21 CFR Part 11 compliance-ready platform

Methods

Sample preparation: The commercially available leuprolide (Figure 1B) was reconstituted in dimethyl sulfoxide. Individual concentrations were prepared by diluting leuprolide in a solvent containing 75:25 (v/v), water/acetonitrile with 1% formic acid and 0.05% human plasma. Individual concentrations were spiked into 100 μ L of human plasma ranging from 2.5 pg/mL to 10000 pg/mL. Protein precipitation was performed with 100 μ L of 5% aqueous ammonia in 90:10 (v/v), acetonitrile/water mixture. Samples were vortexed for 30 seconds and centrifuged at 12000 rcf for 10 minutes at room temperature. The supernatant was transferred to a new Eppendorf tube containing 800 μ L of water and samples were briefly vortexed and centrifuged. Then, 500 μ L of each sample was loaded twice on a Phenomenex Strata-X PRO microelution 96-well plate operated under positive pressure. Samples were washed with 20% acetonitrile in water. Finally, elution was performed using 50 μ L of 1% trifluoroacetic acid in 75:25 (v/v), acetonitrile/water. After elution, samples were diluted with 50 μ L of water after the elution to a final volume of 100 μ L.

Chromatography: Sample separation was performed using an ExionLC system at a 0.6 mL/min flow rate on a Phenomenex Kinetex XB C18 (2.1 x 50 mm, 2.6 μ m, 100 \AA) column. Two cycles of wash in high organic with a flow rate of 1 mL/min were added for column cleaning after the elution of leuprolide (Table 1). A 5-minute gradient was run

using 0.1% formic acid in water as mobile phase A and 0.1% formic acid in acetonitrile as mobile phase B (Table 1). The column temperature was maintained at 55 $^{\circ}$ C. An injection volume of 10 μ L was used for analysis. A mixture of 1:1:1 (v/v/v), acetonitrile/methanol/water was used as a needle wash solvent.

Table 1. Chromatographic gradient for leuprolide.

Time (min)	Mobile phase A (%)	Mobile phase B (%)
0.0	90	10
0.2	90	10
2.5	70	30
2.6*	2	98
2.9*	2	98
3.0	90	10
3.4	90	10
3.5*	2	98
3.9*	2	98
4.0	90	10
5.0	90	10

*Flow rate was set to 1 mL/min

Mass spectrometry: The optimized source and gas parameters are listed in Table 2 and the analyte-dependent MRM parameters are included in Table 3.

Table 2. Source and gas parameters.

Parameter	Value
Polarity	Positive
Ion source gas 1	50 psi
Ion source gas 2	60 psi
Curtain gas	50 psi
Source temperature	500 $^{\circ}$ C
Ion spray voltage	2500 V
CAD gas	12

Table 3. MRM parameters used for quantitation.

ID	Precursor ion (m/z)	Fragment ion (m/z)	CE (V)	CXP (V)	Q0D (V)
Leuprolide	605.6	249.1	40	5	5

Data processing: Data collection and analysis were performed in SCIEX OS software, version 3.0. Peaks were automatically integrated using the MQ4 algorithm and a weighting of $1/x^2$ was used for quantitation.

Quantitative performance

The circulating concentrations in human plasma for depot leuprolide injections vary between high-pg/mL to ng/mL levels.² For lower dose regimens, achieving low-pg/mL sensitivity levels is critical for pharmacokinetic profiling of leuprolide. Here, a sensitive and robust LC-MS/MS quantitation assay with a wide LDR was achieved with a 5-min chromatographic run time.

The calibration curve ranged from 2.5 pg/mL to 10000 pg/mL and was prepared as described in the sample preparation section. Individual concentrations were run in triplicate.

An LLOQ of 2.5 pg/mL was achieved for leuprolide (Figure 1). No interferences were observed in the human plasma matrix blank (Figure 1).

Linearity was achieved between 2.5 pg/mL and 10000 pg/mL with a coefficient of determination (r^2) of 0.993 (Figure 2). An LDR spanning 3.6 orders of magnitude was achieved.

Analytical performance was evaluated for accuracy and precision. The accuracy of the calculated mean was expected to be between 80% and 120% at the LLOQ and between 85% and 115% at higher concentrations. The %CV of the calculated mean for each concentration was expected to be <20% at the LLOQ and <15% at higher concentrations.

Accuracy was within $\pm 12\%$ of the nominal concentration and the %CV was <10% for leuprolide (Figure 3). Calculated accuracy and %CV values met the acceptance criteria at each concentration level.


Figure 2. Calibration curve for leuprolide. The peak area of leuprolide was used to generate a calibration curve. Each concentration level was run in triplicate. Linearity was achieved between 2.5 pg/mL and 10000 pg/mL, spanning an LDR of 3.6 orders of magnitude with an r^2 of 0.993.

Component Name	Actual Concentration	Num. Values	Mean	Standard Deviation	Percent CV	Average Accuracy across Replicates	Value #1	Value #2	Value #3
Leuprolide	2.50	3 of 3	2.650e0	2.620e-1	9.88	106.01	2.502e0	2.496e0	2.953e0
Leuprolide	5.00	3 of 3	4.448e0	7.175e-2	1.61	88.97	4.369e0	4.466e0	4.508e0
Leuprolide	10.00	3 of 3	9.428e0	6.557e-1	6.96	94.28	9.407e0	8.783e0	1.009e1
Leuprolide	25.00	3 of 3	2.602e1	1.625e0	6.24	104.07	2.775e1	2.452e1	2.578e1
Leuprolide	50.00	3 of 3	5.446e1	9.280e-1	1.70	108.93	5.554e1	5.395e1	5.391e1
Leuprolide	100.00	3 of 3	1.018e2	4.910e0	4.82	101.77	9.984e1	9.812e1	1.074e2
Leuprolide	250.00	3 of 3	2.479e2	1.099e0	0.44	99.18	2.492e2	2.472e2	2.474e2
Leuprolide	500.00	3 of 3	5.453e2	2.255e1	4.14	109.05	5.669e2	5.409e2	5.219e2
Leuprolide	1250.00	3 of 3	1.269e3	4.649e1	3.66	101.50	1.321e3	1.252e3	1.233e3
Leuprolide	2500.00	3 of 3	2.476e3	3.136e1	1.27	99.03	2.466e3	2.450e3	2.511e3
Leuprolide	5000.00	3 of 3	4.728e3	9.001e1	1.90	94.57	4.805e3	4.629e3	4.751e3
Leuprolide	10000.00	3 of 3	9.264e3	6.284e1	0.68	92.64	8.192e3	9.304e3	9.298e3

Figure 3. Quantitative performance for leuprolide analysis. Reproducibility and accuracy results were determined from the calibration curve across 3 replicates at each concentration. Statistical results were summarized using the Analytics module in SCIEX OS software.

Compliance-ready SCIEX OS software

SCIEX OS software is a closed system and requires records and signatures to be stored electronically, meeting the regulations outlined by 21 CFR Part 11. SCIEX OS software can open raw data files from any visible storage location within a closed network by using designated processing workstations. Figure 4 illustrates the features of SCIEX OS software that are used for monitoring the audit trail, acquiring and processing data and configuring user access.

The audit trail feature enables users to audit critical user actions and locks in data integrity. The Central Administrator Console (CAC) feature allows users to centralize acquisition and processing using a single platform to maximize efficiency for multi-instrument laboratories, independent of compliance standards. The configuration module allows users to assign roles and access as the administrator, method developer, analyst and reviewer.


Figure 4. Features of the SCIEX OS software for monitoring user access and evaluating the audit trail. The audit trail view allows users to filter for high-risk events easily and enables data integrity features to meet compliance requirements. The software features a Central Administrator Console (CAC) to manage users and groups, role definitions, workstations and projects across all systems. The CAC feature supports both regulated and non-regulated compliance standards. The configuration module enables users to quickly set up roles and levels of access for the administrator, method developer, analyst and reviewer levels.

Conclusions

- An LLOQ of 2.5 pg/mL was reached for the quantitation of leuprolide in human plasma
- Linearity was achieved between 2.5 pg/mL and 10000 pg/mL, with an r^2 of 0.993
- The method demonstrated accurate and highly reproducible (%CV <10%) quantitative performance at all concentrations
- Sensitivity was achieved on the SCIEX 7500 system with an improved front-end technology for better ion generation, capture and transmission
- SCIEX OS software is compliance-ready to support 21 CFR Part 11 and integrates with a triple quadrupole mass spectrometer to support data acquisition, processing and management on a single platform

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An integrated strategy for chemical characterization of Salviaolate lyophilized Injection by SCIEX Triple Quad™ 7500-QTRAP® Ready System

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1. Yaqing Guo, Rui Mao, Yuanyuan Zhang, Rongrong Li, Patrick Kwabena Oduro, Dandan Si, Lifeng Han, Yuhong Huang, Guixiang Pan. An integrated strategy for the systematic chemical characterization of Salvianolate lyophilized injection using four scan modes based on the ultra-high performance liquid chromatography-triple quadrupole-linear ion trap mass spectrometry. *Journal of Pharmaceutical and Biomedical Analysis*, 215 (2022) 114769 doi.org/10.1016/j.jpba.2022.114769;

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RUO-MKT-10-14711-ZH-A

Confident identification of phase 1 metabolites using electron-activated dissociation (EAD)

This technical note describes the identification of pioglitazone phase 1 metabolites from a hepatocyte incubation using comprehensive spectral data from an orthogonal fragmentation technique. More informative MS/MS spectra provided by EAD aided in the software-based identification of phase 1 metabolites to support drug metabolism studies.

Drug metabolism plays a vital role in drug discovery

● Spectrum from UPLC: Pioglitazone EAD 150min, 100°C (pow., 1000) from 3.023 min Precursor: 287.1 Da, +1, CE: 12, C

● Spectrum from MetS: Pioglitazone CID 150min, 100°C (pow., 1000) from 3.023 min Precursor: 287.1 Da, +1, CE: 15, C


Figure 1. EAD provided rich MS/MS spectra for confident identification of pioglitazone phase 1 metabolites. EAD (top) and CID (bottom) spectra for a metabolite showing aliphatic hydroxylation and alcoholic oxidation at a retention time of 3.77 minutes.

and development, affecting pharmacokinetics, pharmacodynamics and safety. Studies of in vitro metabolism of drugs in human and animal tissues help identify major metabolism pathways ("soft spots").¹

The Mass-MetaSite software employs several algorithms to detect the peaks corresponding to metabolites in the mass spectra obtained from incubated samples. This process includes background subtraction, noise suppression, isotope pattern analysis, retention time analysis and mass shift analysis based on cyp, non-cyp and other uncommon reactions. A potential metabolite is scored based on the number of matches between its fragments and the parent compound. Site of metabolism (SoM) predictions generated by the Mass-MetaSite software can be used to distinguish between potential regioisometric metabolites that have the same fragmentation patterns and mass shifts. This technical note demonstrates an efficient soft-spot identification workflow using a novel orthogonal fragmentation technique, EAD, on the ZenoTOF 7600 system (Figure 1). Sites of metabolism were predicted using the Mass-MetaSite software.

Key features for the identification of metabolites using the ZenoTOF 7600 system

- **Confident identification:** Acquire more fragments to confidently identify the site of metabolism for phase 1 metabolites with EAD. EAD spectra were more informative than collision-induced dissociation (CID) spectra for metabolite identification.
- **Fast characterization and identification:** Achieve rapid and efficient software-aided characterization and identification of drug metabolites from hepatocyte

incubations using the ZenoTOF 7600 system

- **Detection of low-level metabolites:** Identify low-level metabolites present in drug metabolism studies with enhanced MS/MS sensitivity provided by the Zeno trap
- **Streamlined data acquisition and processing workflow:** Develop confident structure-metabolic stability relationships for drugs using a quick, easy-to-use workflow from acquisition to analysis

Methods

Sample preparation: Pioglitazone was incubated at 37 °C in human hepatocytes at a starting concentration of 5µM. Samples were removed from incubation and quenched with acetonitrile at 0-, 30-, 60-, 90- and 120-minute intervals.

Chromatography: Separation was performed on a Phenomenex Kinetex Polar C18 column (2.1 x 100 mm, 2.6 µm, 100 Å) at a column temperature of 40°C. Mobile phase A was 0.1% (v/v) formic acid in water and mobile phase B was 0.1% (v/v) formic acid in methanol. An injection of 5 µL was subjected for analysis.

The chromatographic gradient conditions used are summarized in Table 1.

Table 1. Chromatographic gradient.

Time (min)	Mobile phase A (%)	Mobile phase B (%)
0.0	95	5
1.0	95	5
7.0	5	95
9.0	5	95
9.1	95	5
10	95	5

Mass spectrometry: The samples were analyzed using the data-dependent acquisition (DDA) method with Zeno CID DDA and Zeno EAD DDA on the ZenoTOF 7600 system. The source and gas conditions used are summarized in Table 2.

The method conditions are summarized in Table 3.

Table 2. Source and gas conditions.

Parameter	Setting
Curtain gas	40 psi
Ion source gas 1	55 psi
Ion source gas 2	65 psi
CAD gas	7
Ion spray voltage	5500 V
Source temperature	500 °C

Table 3. Zeno DDA parameters.

Parameter	Setting
Method duration	10 min
TOF MS start-stop mass	100–1000 Da
Maximum candidate ions	5
Accumulation time (TOF MS)	0.1 s
TOF MS/MS start-stop mass	50–1000 Da
Accumulation time (TOF MS/MS)	0.1 s
Collision energy (CID)	40 V
Collision energy spread (CID)	15 V
Electron kinetic energy (EAD)	12 eV
Electron beam current (EAD)	6000 nA

Data processing: SCIEX OS software, version 3.0 was used for data acquisition. The Mass-MetaSite software was used to predict biotransformation sites using Zeno CID DDA and Zeno EAD DDA data.⁴⁻⁹

EAD provides positional information to identify the site of metabolism

In the 120-minute incubation sample, 3 peaks at retention times of 3.53, 3.77 and 4.13 minutes were identified as hydroxy pioglitazone.

The Mass-MetaSite software identified the peaks as hydroxy pioglitazone using both Zeno CID DDA and Zeno EAD DDA

data. The Zeno EAD data for the hydroxy pioglitazone peak at 3.53 minutes yielded 10 product ion matches. Therefore, the software labeled this peak hydroxy pioglitazone (M-VII). In comparison, for the Zeno CID data, the software predicted 2 possibilities for hydroxylation, which included hydroxy pioglitazone (M-II) and hydroxy pioglitazone (M-IV). These 2 structures were predicted with equal likelihoods and 7 product ion matches each (Figure 2).


Figure 2. EAD and CID spectra for the hydroxy pioglitazone metabolite at retention time 3.53 minutes with fragment ion matches for structures predicted by Mass-MetaSite software. Product ion matches with pioglitazone are displayed in red and metabolite-specific matches are indicated in yellow.

The hydroxy pioglitazone peak at retention time 3.77 minutes was scored highest for hydroxy pioglitazone (M-IV) based on 6 product ion matches. The hydroxy pioglitazone (M-II) was ranked second with 6 fragment matches from

CID data. EAD data resulted in 6 possible structures for the hydroxyl metabolite. The hydroxy pioglitazone (M-IV) was ranked first with the highest score and 7 product ion matches (Figure 3).


Figure 3. EAD spectra for the hydroxy pioglitazone metabolite at retention time 3.77 minutes with fragment ion matches for structures predicted by Mass-MetaSite software and SoM prioritization ranking. Product ion matches with pioglitazone are displayed in red and metabolite-specific matches are indicated in yellow. The SoM prioritization ranking indicated the highest probability match for hydroxy pioglitazone (M-IV).

The hydroxy pioglitazone peak at retention time 4.13 minutes was labeled as a possible amide hydrolysis and dehydrogenation based on 6 product ion matches or hydroxy pioglitazone (M-II) with low probability based on 5 product ion matches from the CID data. The EAD data

indicated a higher likelihood that the peak corresponded to hydroxy pioglitazone (M-II) based on 8 product ion matches (Figure 4).


Figure 4. EAD and CID spectra for the hydroxy pioglitazone metabolite at retention time 4.13 minutes with fragment ion matches for structures predicted by Mass-MetaSite software. Product ion matches with pioglitazone are displayed in red and metabolite-specific matches are indicated in yellow.

A metabolite with aliphatic hydroxylation and alcoholic oxidation was labeled at retention time 3.78 minutes using both EAD and CID spectra. The rich MS/MS spectra collected

by EAD enabled a more confident structure assignment. The EAD data yielded 8 fragment matches, whereas the CID data yielded only 3 fragment matches (Figure 5).³⁻⁹


Figure 5. EAD and CID spectra for a metabolite showing aliphatic hydroxylation and alcoholic oxidation with fragment ion matches for structures predicted by Mass-MetaSite software. Product ion matches with pioglitazone are displayed in red and metabolite-specific matches are indicated in yellow.

All 3 pioglitazone hydroxy metabolites (M-VII, M-IV and M-II) were confirmed by matching their respective retention time with standard solution injections. The ZenoTOF 7600 system demonstrated excellent mass accuracy for the workflow. All metabolites and fragments were identified with <2 ppm error, which enabled the confident identification of phase 1 metabolites present in an in vitro metabolism study of pioglitazone. Furthermore, the identification of phase 1 metabolites was easily and confidently performed with the MS/MS coverage provided by Zeno EAD on the ZenoTOF 7600 system.

Conclusions

- More information-rich product ion spectra provided by EAD aided in the rapid software-aided characterization and identification of phase 1 metabolites from hepatocyte incubations of pioglitazone on the ZenoTOF 7600 system
- The presented workflow can be easily adapted for in vivo metabolism studies for the detection of low-level metabolites with the enhanced sensitivity provided by the Zeno trap
- A streamlined data acquisition and processing workflow was utilized to expedite data reduction and develop confident structure-metabolic stability relationships for pharmaceutical drugs

References

1. Zhang, Z.; Tang, W. (2018). Drug metabolism in drug discovery and development. *Acta Pharm Sin B*. 8(5): 721-732. Orthogonal fragmentation mechanism enables new levels of metabolite characterization, SCIEX technical note, RUO-MKT- 02-13348-A.
2. Boon, B. et al. (2010). Enhanced metabolite identification with MS(E) and a semi-automated software for structural elucidation. *Rapid Commun Mass Spectrom*. 24(21): 3127-38.
3. Zamora, I. et al. (2013). High-throughput, computer-assisted, specific MetID. A revolution for drug discovery. *Drug Discovery Today Technol*. 10(1): e199-205.
4. Strano-Rossi, S. et al (2014). Metabolism of JWH-015, JWH-098, JWH-251, and JWH-307 in silico and in vitro: a pilot study for the detection of unknown synthetic cannabinoids metabolites. *Analytical and Bioanalytical Chemistry*. 406(15): 3621-3636.
5. Zelesky, V. (2013). Software automation tools for increased throughput metabolic soft-spot identification in early drug discovery. *Bioanalysis*, 5(10): 1165-1179.
6. Li, A.C. (2013). Update on hydrocodone metabolites in rats and dogs aided with a semi-automatic software for metabolite identification *Mass-MetaSite*. *Xenobiotica*, 43(4): 390-398.
7. Brink, A. et al. (2014). Post-acquisition analysis of untargeted accurate mass quadrupole time-of-flight MS(E) data for multiple collision-induced neutral losses and fragment ions of glutathione conjugates. *Rapid Commun Mass Spectrom*. 28(24): 2695-703.
8. Ge, S; Tu, Y; Hu, M. Challenges and opportunities with predicting in vivo phase II metabolism via glucuronidation from in vitro data (2016). *Curr Pharmacol Rep*. 2(6): 326-338

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Using The EAD Technology of SCIEX ZenoTOF® 7600 system to Identify leucine and Isoleucine in TESTUDINIS CARAPACIS ET PLASTRI COLLA

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Xiaoyuan Shi, Chenchun Zhong, Zhimin Long, Lihai Guo

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© 2023 DH Tech. Dev. Pte. Ltd. MKT-29782-ZH-A